



Operating Instructions

WHEELCHAIR ACCESS VEHICLE WITH REAR CUT-OUT

RENAULT KANGOO / MERCEDES-BENZ CITAN WITH AMF-BRUNS EASY-RAMP AND EASY-FLEX-RAMP



www.amf-bruns.de

Foreword

Dear reader,

these Operating Instructions provide all information necessary to operate the Easy and Easy-Flex wheelchair access ramps safely.

The access ramp is designed and constructed in accordance with state of the art technology and recognised safety standards. Persons and materials can however be endangered, as not all danger areas can be eliminated if the functional capability is to be maintained. Accidents caused by these dangers can however be prevented by strictly observing these Operating Instructions. Over and above this, the operational efficiency of the access ramp can be used to the full and unnecessary faults can be prevented.

After reading these Operating Instructions for the first time, keep them in a safe place for future reference over the entire lifetime of the vehicle. If you sell the vehicle, hand these Operating Instructions over to the new owner.

All details, figures and dimensions contained in these Operating Instructions are non-binding. No claims in any form can be derived from these.

This document must not be reproduced or duplicated, in full or in part, without the prior, written permission of the manufacturer.

The access ramp must never be converted or modified in any way, without seeking the prior, written permission of the manufacturer. The manufacturer will not be held responsible in any way whatsoever if conversions or modifications are carried out without authorisation.

Use only original spare parts or spare parts which have been approved of by the manufacturer. If spare parts other than these are used, this can have a negative effect on the specified characteristics, the functionality and safety of the access ramp. Using non-original or unauthorised spare parts will render the guarantee null and void.

Contact our customer services department to order spare parts or accessories (see Chapter 10, page 47).

Explanation of symbols and signs

To improve understanding, the following conventions should be met for these Operating Instructions:

1.

The following types of special note are used to highlight text passages of an important nature.



DANGER!

- warns of a situation of immediate danger, which will lead to serious or fatal injuries, if not avoided.



WARNING!

- warns of a potentially dangerous situation, which will lead to serious or fatal injuries, if not avoided.



CAUTION!

- warns of a potentially dangerous situation, which will lead to slight or minor injuries, or material damage if not avoided.



ATTENTION!

...warns of a potentially dangerous situation, which can cause material damage, if not avoided.



...contains general notes and useful information.



...gives a reference to important information in other sections and documents.

2.

Some texts serve a particular purpose. These are identified as follows:

- Lists

⇒ Instructional text, e.g. a sequence of activities.

3.

Meaning of directions:

If directions are given in the text (in front of, front, behind, rear, right, left), these directions relate to the normal direction of travel of the vehicle.



Contents

| | |
|--|-----------|
| Foreword..... | 3 |
| Contents..... | 6 |
| 1 Safety | 8 |
| 1.1 Proper Use..... | 8 |
| 1.2 Improper Use | 9 |
| 1.3 Personnel Requirements..... | 10 |
| 1.4 Product Monitoring | 10 |
| 1.5 Danger Zone | 11 |
| 1.6 Safety Devices | 11 |
| 1.7 Safety and Accident Prevention Regulations | 12 |
| 1.8 De-Commissioning and Conservation | 13 |
| 1.9 Disposal | 13 |
| 2 Description | 14 |
| 2.1 Layout and Function..... | 14 |
| 2.1.1 Floor section..... | 15 |
| 2.1.2 Ramp..... | 15 |
| 2.1.3 Belts and belt anchorages..... | 16 |
| 2.2 Rating Plate | 17 |
| 2.3 Operating Controls..... | 18 |
| 2.4 Technical Data | 19 |
| 3 Transportation..... | 19 |
| 4 Installation / Commissioning | 19 |
| 5 Operation..... | 20 |
| 5.1 Safety regulations for operation | 20 |
| 5.2 Deploying the Ramp..... | 22 |
| 5.3 Loading a Passenger / Wheelchair..... | 24 |
| 5.4 Disembarking a Passenger / Unloading a Wheelchair | 27 |
| 5.5 Easy-Flex: Stowing the Ramp | 29 |

| | |
|--|-----------|
| 6 Maintenance and Repair | 31 |
| 6.1 Safety Regulations for Maintenance and Repair | 31 |
| 6.2 Routine Maintenance Work | 32 |
| 6.2.1 Maintenance schedule | 32 |
| 6.2.2 Maintenance record | 32 |
| 6.3 Functionally Testing the Electrically Operated Belts..... | 33 |
| 6.4 Yearly Inspection..... | 34 |
| 6.5 Maintenance and Inspection Record..... | 34 |
| 7 Faults and Troubleshooting | 35 |
| 7.1 Troubleshooting Table..... | 35 |
| 8 Inspection Log Book | 38 |
| 8.1 Inspection Log Book Master Data Sheet..... | 39 |
| 8.2 Inspection List..... | 40 |
| 8.3 Inspection Results..... | 41 |
| 9 Electrical Circuit Diagram | 46 |
| 10 Customer Service | 47 |

1 Safety



CAUTION!

There are a number of risks of suffering personal injury and material damage involved in the operation and maintenance of the access ramp.

Therefore:

- it is imperative, that these Operating Instructions are read thoroughly before operating your access ramp. Always observe the notes and information contained therein, in particular the safety instructions.
 - If these Operating Instructions or parts thereof are lost or become illegible, please request a new copy from the manufacturer.
-

Prerequisite to the safe handling and trouble-free operation of the access ramp is a thorough knowledge of the applicable safety information and the safety regulations.

It is therefore imperative that this Chapter is read thoroughly before operating the access ramp and that the instructions and warnings herein are strictly observed. The safety information and warnings, given at the appropriate places in the following Chapters, must also be strictly observed. The manufacturer will not be held responsible if safety information and warnings are not observed.

In addition to the information given in these Operating Instructions, local legislative regulations must be taken into consideration, in particular those regarding safety and accident prevention.

1.1 Proper Use

The access ramp must only be used for transporting a disabled person in a wheelchair or an empty wheelchair. By using the access ramp, persons can be loaded into or disembarked from the vehicle to which it is fitted.

When doing so, the access ramp will be operated by an accompanying person.

Wheelchairs must be equipped with the fittings necessary to secure it by means of the belt restraint systems provided.

Proper use also includes strictly adhering to the information given in these Operating Instructions.



WARNING!

If the access ramp is used for any other purpose than that described above, this may result in dangerous situations for persons or material damage being caused.

Therefore:

- only use the access ramp for the purpose for which it was intended.
 - Always adhere to information given in these Operating Instructions.
 - Do not use the access ramp for any other use, particularly those given in Section 1.2. These are deemed to be improper use.
-

1.2 Improper Use

Any type of use, other than that mentioned in Section 1.1 is deemed to be improper use.

The access ramp is deemed to be improperly used if for example:

- it is used to transport goods,
- it is used in a faulty state or with safety-relevant malfunctions,
- it is used to transport wheelchairs that can not be secured by means of the belt restraint systems provided, or
- it is used by incompetent persons.

1.3 Personnel Requirements

The access ramp must only be handled by persons who:

- have been instructed in how to operate the access ramp,
- have read and understood these Operating Instructions,
- have proven their ability to operate the access ramp to the vehicle's owner,
- have been expressly assigned by the vehicle's owner to operate the access ramp,
- have the technical knowledge to operate the wheelchair's brakes and switch the motors of electrically driven wheelchairs ON and OFF and
- are in a position to adapt themselves to the particular behaviour and needs of disabled persons.

Transportation, installation, commissioning, maintenance, repair, fault finding and disposal of the access ramp and the vehicle itself must only be carried out by persons with the corresponding technical training and experience.

1.4 Product Monitoring

Please contact AMF-Brunns GmbH & Co. KG immediately if faults or problems are encountered when operating your access ramp or if accidents or "near-misses" occur.

AMF-Brunns will effect a solution to the problem with your help and the knowledge gained will flow into future projects.



NOTE

Guarantee work on the access ramp must only be carried out with the prior agreement of AMF-Brunns GmbH & Co. KG.

The costs of such work will not be accepted by AMF-Brunns without prior agreement.

1.5 Danger Zone

The danger zone is the area behind the access ramp that is being deployed, in which a person can be injured through the movement of the access ramp.

CAUTION!

Risk of injury and danger of material damage when deploying the access ramp.

There is a risk of injuries to the legs and feet if the access ramp is allowed to fall to the ground out of control. The access ramp can be damaged.

Therefore:

- keep hold of the access ramp until it is fully in contact with the ground.
 - Allow the access ramp to deploy slowly.
-



1.6 Safety Devices

Electrically operated belts

Electrically operated belts, with which the wheelchair is secured when driving, are attached to the front area of the vehicle's floor. These electrically operated belts can only be pulled out when the ramp has been fully deployed and the indicator lamp in the rocker switch is lit (refer to Section 2.3, page 18). The electrically operated belts retract automatically when the rocker switch is in either position. Before the wheelchair can be pushed or driven up the ramp, the electrically operated belts must be attached to the wheelchair and the rocker switch must be operated such that the indicator lamp goes out. This will prevent the wheelchair from rolling back out of control.

1.7 Safety and Accident Prevention Regulations

Adhere to the following notes in order to prevent personal injuries and material damage. Adhere also to the relevant safety and accident prevention regulations laid down by the trade associations for commercial use.

- The access ramp must only be operated if all safety devices are correctly fitted (see Section 1.6, page 11). Such devices must only be removed in order to enable maintenance and repair work to be carried out. All safety devices must be replaced immediately after such work has been completed. If they are not replaced correctly, there is a high risk of injury.
- The access ramp must only be used for the purpose for which it is intended, otherwise dangerous situations, with resultant injuries, may occur (Proper use: see Section 1.1, page 8).
- The owner is responsible for ensuring that proper use is adhered to, in particular that the access ramp is only operated by authorised persons.
- If the access ramp is used commercially or as a public utility, the owner must ensure that operating personnel are familiar with the operation of the access ramp under all operating conditions by giving training and familiarisation courses.
- Proper use of the access ramp also includes adherence to the specified maintenance and repair work, in particular strict adherence to the maintenance intervals. If such work is not carried out, trouble-free operation can not be guaranteed. Risks of personal injury and material damage can occur. We recommend that maintenance records be kept.
- If the access ramp is used commercially or as a public utility, it must be inspected by a technical expert at intervals of not more than 1 year after commissioning. During inspection, faults affecting the safety should be systematically identified and remedial action taken (see "Inspection Log Book", page 36).

- The access ramp must not be operated in a faulty condition, as severe injuries can be caused by this. If faults occur, do not use the access ramp until repairs have been effected.
- Before carrying out maintenance or repair work, this includes cleaning work, switch the vehicle's engine OFF and prevent the vehicle from rolling away by applying the handbrake. Make certain that no other person can start the vehicle (e.g. by removing the ignition key). If this is not done, there is a risk of injury.
- Use only original spare parts or spare parts which have been approved of by the manufacturer. If non-original or unauthorised replacement parts are used, this will render the guarantee null and void.

1.8 De-Commissioning and Conservation

For queries regarding de-commissioning and conservation, contact our customer services department (see Chapter 10, page 47).

1.9 Disposal

When the access ramp and vehicle's useful life has expired, they must only be disposed of by qualified specialists. The manufacturer will accept no liability for damage caused by incorrect disposal.

2 Description

Persons who are wheelchair-bound can be loaded into or disembarked from a vehicle via the access ramp.

The main components of the access ramp are:

- the floor section,
- the ramp,
- the electrically operated belts and retractors for securing the wheelchair.

The aim of this chapter is to illustrate the construction and function of the access ramp. To this end, the individual assemblies and components are described in the following sections.

2.1 Layout and Function



Figure 1: Ramp in The Deployed Position

2.1.1 Floor section

The floor section (see Figure 1, page 14) offers sufficient space for the passenger in the wheelchair. The surface of the floor section is lower than the floor of the vehicle. This ensures that the passenger has sufficient headroom.

2.1.2 Ramp

The Easy-Ramp is used in two positions.

In the completely deployed position, the ramp provides a slope, via which the wheelchair can be easily pushed or driven into the vehicle (see Figure 1, page 14).

The rear doors can be closed when it is in the vertical position. The ramp leaves sufficient space for a passenger or an empty wheelchair. In this position, the frame of the access ramp fills the gap in the vehicle's rear bumper (see Figure 2).




Figure 2: Easy-Ramp in the Vertical Position

A third position is available for the Easy-Flex-Ramp. The ramp can be folded down into the vehicle, so that it is at the same level as the floor of the vehicle (see Figure 3, page 16).

The floor section is almost completely covered. This creates a level luggage compartment floor. The vehicle can be used without difficulty to transport goods and other items, up to a weight of 100 kg, in the normal way.



Figure 3: Fully Stowed Easy-Flex-Ramp

2.1.3 Belts and belt anchorages

Two electrically operated belts are fitted to the front area of the floor section (see Figure 4). The electrically operated belts can be pulled out when the ramp is fully deployed and the indicator lamp in the rocker switch is lit.




Figure 4: Electrically operated belts

Fittings, to which the retractor belts can be attached, are provided on the rear area of the floor section (see Figure 5). The rear end of the wheelchair is secured by the retractor belts.




Figure 5: Retractor Belts and Fittings

2.2 Rating Plate

A rating plate, which contains the fundamental data, is attached to the access ramp (see Figure 6).




Figure 6: Rating Plate

2.3 Operating Controls




Figure 7: Rocker Switch




Figure 8: Operating Levers

| Designation | Function |
|--|--|
| Rocker switch | Allows the electrically operated belts to be pulled out |
| Operating lever 1 (only for Easy-Flex) | Allows the ramp to be folded down inside the vehicle to the horizontal position. |
| Operating lever 2 | Allows the ramp to be deployed. |

2.4 Technical Data

| | |
|--|-----------------------------|
| Designation | Easy-Ramp Easy-Flex-Ramp |
| Maximum load (wheelchair use) | 300 kg |
| Maximum load (Easy-Flex-Ramp stowed position, luggage compartment) | 100 kg |

3 Transportation

Transportation of the access ramp is carried out by the manufacturer or by trained, experienced personnel.

4 Installation / Commissioning



The access ramp must be installed in accordance with the AMF-Brunns Installation Instructions applicable to the vehicle in question.

DANGER!

Risk of injury through incorrect installation.

A number of risks of personal injury and material damage can be caused if the access ramp is incorrectly installed in the vehicle.

Such risks of danger can not only occur during installation but also as a result of installation not being carried out correctly.

Therefore:



- the access ramp must only be installed by specialist personnel, who have been trained for this job by the manufacturer. This also applies to the electrical connections.
- The vehicle manufacturer's body fitting guidelines must be adhered to.
- A technical expert must be called in to inspect the access ramp for correct installation and the effectiveness of the safety devices.
- The technical expert must confirm the fact that operational safety of the access ramp has been established by making a corresponding entry in the inspection log book (see Chapter 7.2, page 36).
- The access ramp must not be used until this has been done.

5 Operation

5.1 Safety regulations for operation



Before operating the access ramp, the safety information must have been read (see Chapter 1, page 8).



WARNING!

Risk of injury and material damage if the access ramp is operated in a faulty condition.

Therefore:

- Always carry out inspections in accordance with the maintenance schedule before use (see Chapter 6, page 31).
- Have the access ramp inspected at intervals of not more than one year by a technical expert.



WARNING!

Risk of injury and material damage if the access ramp is operated by unauthorised persons.

Dangerous operating conditions may be caused if the access ramp is operated by unauthorised persons.

Therefore:

- The access ramp must only be operated by persons who are familiar with operating the access ramp.
- The access ramp must never be operated by the passenger.
- Lock the vehicle's doors when the access ramp is not in use.



WARNING!

Risk of injury if safety devices are removed or are defective.

If safety devices are modified, bypassed or removed, they will no longer fulfil their function.

Therefore:

- Never modify, bypass or remove safety devices.
- Always refit safety devices if they have been removed (e. g. for maintenance or repair purposes).



WARNING!

Risks through road traffic collisions.

A number of risks are involved for the passenger operator when using the access ramp on public highways.

Therefore:

- When parking the vehicle ensure that other road traffic is affected as little as possible.
 - Park the vehicle such that as little danger as possible is presented to passenger and operator by other road traffic.
 - Switch the vehicle's engine OFF and apply the handbrake before operating the access ramp.
-



CAUTION!

Risk of injury if there is insufficient space to operate the access ramp.

If there is insufficient space to leave or to get onto the access ramp, there is a risk of being injured by the access ramp or other objects (e. g. walls, posts or poles).

These risks apply to both passenger and accompanying person.

Therefore:

- When parking the vehicle, make certain that there is sufficient space to operate the access ramp.
-



CAUTION!

Risk of personal injury and material damage.

The rear parking distance control warning characteristic is changed by fitting the ramp. The parking distance control first warns at a shorter distance from obstacles than usual.

Therefore:

- Bear the change in warning characteristic in mind when reversing.
-

5.2 Deploying the Ramp

- ⇒ Switch the vehicle's engine OFF.
- ⇒ Remove the ignition key.
- ⇒ Apply the vehicle's handbrake.
- ⇒ Open the rear doors.
- ⇒ For the Easy-Flex-Ramp: Use the handhold to fold the ramp back up to the vertical position, until it audibly engages (see Figure 9).

CAUTION!

Risk of crushing when lifting the ramp to the vertical position. When lifting the ramp to the vertical position, there is a risk of being crushed between the ramp and the frame (rear bumper).

Therefore:

- Do not reach into the gap between ramp and frame.



Figure 9: Folding the Ramp to the Vertical Position (only Easy-Flex)

- ⇒ Check that the ramp is securely locked in the vertical position by gently shaking it.
- ⇒ Keep a firm hold on the ramp.
- ⇒ Push operating lever 2 to the left (see Figure 10, page 23).



Figure 10: Unlocking the Frame

⇒ Fold the ramp out until it makes contact with the ground (see Figure 11, page 24).

CAUTION!

Risk of injury and material damage when deploying the ramp.

There is a risk of injuries to the legs and feet if the ramp is allowed to fall to the ground out of control. The ramp can also be damaged.

Therefore:

- Keep a firm hold of the ramp until it is fully in contact with the ground.
- Allow the ramp to deploy slowly.





Figure 11: Deploying the Ramp

5.3 Loading a Passenger / Wheelchair

- ⇒ Deploy the ramp (see Section 5.2, page 22).
- ⇒ Accompany the wheelchair to a position in front of the deployed ramp.
- ⇒ Operate the rocker switch such that the indicator lamp in the rocker switch lights up (see Figure 7, page 18).
- ⇒ Pull the electrically operated belts out to the wheelchair.
- ⇒ Attach the electrically operated belts to the wheelchair (see Figure 12, page 25).



Figure 12: Attaching the Electrically Operated Belts to the Wheelchair

⇒ Operate the rocker switch such that the indicator lamp in the rocker switch goes out.

WARNING!

Risk of injury through the wheelchair rolling backwards.

The floor section and the ramp form a slope via which the wheelchair can roll back out of control. The passenger in the wheelchair can be severely injured. Material damage can be caused.



Therefore:

- Attach the electrically operated belts to the wheelchair.
- Operate the rocker switch such that the indicator lamp in the rocker switch goes out.
- Push or drive the wheelchair up the access ramp only when this has been done.

⇒ Push or drive the wheelchair into the vehicle.

WARNING!

Risk of injury if the electrically operated belts loosen when driving.

If the electrically operated belts are unevenly tensioned, they can become loose when driving. There is a risk of the passenger in the wheelchair and of the other occupants in the vehicle being injured



Therefore:

- Push the wheelchair into the centre of the vehicle.
- Make sure that the electrically operated belts leave the belt reels at the same angle as one another.
- If this is not the case, re-position the wheelchair.

- ⇒ Switch the motor of an electrically driven wheelchair OFF.
- ⇒ Attach the retractor belts to the wheelchair.
- ⇒ Clip the retractor belts onto the fittings in the rear area of the floor of the vehicle (see Figure 5, page 17).
- ⇒ Tension the retractor belts by tightening the belt reels (see Figure 13).




Figure 13: Tensioning the Retractor Belts

- ⇒ Apply the wheelchair's brake(s).
- ⇒ Secure the passenger in the wheelchair with the lap belt, the automatic diagonal shoulder belt and the head and back support FutureSafe (optional).

⇒ Fold the ramp back up to the vertical position, until it audibly engages (see Figure 14).

CAUTION!

Danger of crushing when folding the ramp to the vertical position.



When lifting the ramp to the vertical position, there is a risk of hands or fingers being crushed between the ramp and the frame (rear bumper).

Therefore:

- Do not reach into the gap between the ramp and the frame.



Figure 14: Folding the Ramp Up

⇒ Close the rear doors.

5.4 Disembarking a Passenger / Unloading a Wheelchair

- ⇒ Deploy the ramp (see Section 5.2, page 22).
- ⇒ Swivel the FutureSafe head and back support (optional) out of the vehicle.
- ⇒ Release the lap belt and the automatic diagonal shoulder belt with which the passenger is secured (optional).
- ⇒ Slacken the retractor belt reels (see Figure 13, page 26).

- ⇒ Unclip the retractor belts from the fittings (see Figure 5, page 17).
 - ⇒ Detach the retractor belts from the wheelchair.
 - ⇒ Switch the motor of an electrically driven wheelchair ON.
 - ⇒ Release the wheelchair's brake(s).
 - ⇒ Take up a position behind the wheelchair so that you can stop the wheelchair safely if it rolls backwards.
-

WARNING!

Risk of injury through a wheelchair rolling backwards out of control.



The floor section and the ramp form a slope via which the wheelchair can roll down out of control. The passenger in the wheelchair can be severely injured. Material damage can be caused.

Therefore:

- Stand behind the wheelchair to prevent it from rolling backwards.
 - Only then operate the rocker switch such that the indicator lamp in the rocker switch lights up.
-

- ⇒ Operate the rocker switch such that the indicator lamp in the rocker switch lights up (see Figure 7, page 18).
- ⇒ Accompany the passenger out of the vehicle.
- ⇒ Detach the electrically operated belts from the wheelchair.
- ⇒ Allow the electrically operated belts to retract.
- ⇒ Operate the rocker switch such that the indicator lamp in the rocker switch goes out (see Figure 7, page 18).
- ⇒ Fold the ramp back up to the vertical position, until it audibly engages (see Figure 14, page 27).
- ⇒ Close the rear doors.

5.5 Easy-Flex: Stowing the Ramp

With the access ramp fully folded away and the rear seats in their original position, the vehicle can be used as a standard vehicle (see Figure 15).



Figure 15: Ramp Fully Folded Down (only Easy-Flex)

⇒ Push operating lever 1 to the right (see Figure 16 and Figure 8, page 18).



Figure 16: Unlatching the Ramp



CAUTION!

Risk of injury and danger of material damage when stowing the ramp fully in the vehicle.

When folding the ramp to the horizontal position, there is a risk of being crushed between the access ramp and the floor section. The access ramp can be damaged.

Therefore:

- Keep a firm hold of the ramp until it is fully in contact with the floor section.
- Do not reach into the gap between access ramp and floor section.

⇒ Fold the ramp down slowly, until it's bearing surfaces rest on the vehicle's floor at the sides of the floor section (see Figure 17).



Figure 17: Stowing the Ramp

The ramp rests on both sides of the floor section and can be used as a normal luggage compartment floor (see Figure 15, page 29).



ATTENTION!

The ramp can suffer damage if it is overloaded.

Therefore:

- Do not stand on the stowed access ramp.
- Do not place a load of more than 100 kg on the stowed access ramp.

6 Maintenance and Repair

6.1 Safety Regulations for Maintenance and Repair



Read the safety information before carrying out maintenance and repair work (see Chapter 1, page 8).



WARNING!

Risk of injury through missing safety devices.

If safety devices are removed, they will no longer fulfil their function.

Therefore:

- Always refit safety devices if they have been removed (e. g. for maintenance or repair purposes).
-



WARNING!

Risk of injury and material damage if maintenance and repair work is not carried out correctly.

Therefore:

- Only allow specialist personnel to carry out maintenance and repair work.
 - If repair work has been carried out on load bearing parts, this must be inspected by a technical expert.
-



WARNING!

Risk of personal injury or material damage if inferior quality spare parts are used.

The manufacturer will accept no liability whatsoever for damage or injury caused by the use of non-original spare parts or spare parts that have not been approved of by the manufacturer.

Therefore:

- Use only original spare parts or spare parts that have been approved of by the manufacturer.
-

6.2 Routine Maintenance Work

The maintenance schedule below lists the maintenance work that must be carried out at regular intervals.

Contact the customer service department for maintenance work that must be carried out by a specialist (see Chapter 10, page 47).

6.2.1 Maintenance schedule

| Interval | Activity | see |
|-------------|--|-------------|
| Before use | Functionally test the electrically operated belts. | Section 6.3 |
| | Carry out a visual inspection for damage or faults. Initiate repairs if necessary. | |
| yearly | Have the access ramp inspected by a technical expert. | Section |
| as required | Clean the access ramp. | |

6.2.2 Maintenance record

Enter maintenance and repair work that has been carried out into the maintenance record provided for this purpose (see Section 6.5, page 34). This provides a traceable record of maintenance work.

For records of maintenance work over and above this, it is recommended that you keep your own lists.

6.3 Functionally Testing the Electrically Operated Belts

- ⇒ Fold the access ramp up to the vertical position, until the spring-loaded locking bolts audibly engage.
- ⇒ Operate the rocker switch several times (see Figure 7, page 18). Observe whether the indicator lamp in the red rocker switch lights up when doing so.
- ⇒ Fold the access ramp out until it rests on the ground.

CAUTION!

Risk of injury and danger of material damage when deploying the ramp.



There is a risk of injuries to the legs and feet if the ramp is allowed to fall to the ground out of control. The ramp can be damaged.

Therefore:

- Keep hold of the access ramp until it is fully in contact with the ground.
- Allow the ramp to deploy slowly.

-
- ⇒ Operate the rocker switch such that the indicator lamp in the rocker switch goes out.
 - ⇒ Try to pull the electrically operated belts out.
 - ⇒ Operate the rocker switch again such that the indicator lamp lights up.
 - ⇒ Try again to pull the electrically operated belts out.

The indicator lamp in the rocker switch must only light up when the ramp is deployed. The electrically operated belts must only be able to be pulled out when the access ramp is deployed and the indicator lamp in the rocker switch is lit.

If this is not the case:

- ⇒ Do not use the access ramp and have the access ramp's electrical system checked.

7 Faults and Troubleshooting



WARNING!

Risk of injury and material damage if repair work is carried out incorrectly.

Therefore:

- Only allow specialist personnel to carry out repair work.

If faults occur when operating the access ramp, proceed as described in the following troubleshooting table. Contact the customer service department if faults are encountered which cannot be remedied using the information and measures given in the table.

7.1 Troubleshooting Table

| Fault | (possible) Cause | Remedial Measures |
|--|---|---|
| The ramp does not deploy or stow correctly. | Foreign objects blocking the latching mechanisms. | Inspect the latching mechanisms and remove foreign objects, clean the latching mechanisms if necessary. |
| | A spring-loaded latch is defective. | Check the spring-loaded latches and replace if necessary. |
| | One or more of the gas-operated springs is/are defective. | Check the gas-operated springs and replace them if necessary. |
| | The access ramp has too little play. | Adjust the access ramp play (see Section 7.2, page 36) |
| The access ramp rattles when driving. | The access ramp has too much play. | Adjust the access ramp play (see Section 7.2, page 36) |
| The electrically operated belts can not be pulled out. | The rocker switch is defective. | Have a new rocker switch fitted. |
| | The fuse has blown. | Replace the fuse in the vicinity of the battery. |
| | The electrically operated belts are defective. | Have the electrically operated belts replaced. |

7.2 Adjusting the access ramp's play

The access ramp can show signs of play when in the vertical position. This causes the access ramp to rattle when driving. Adjusting screws are provided at both sides of the access ramp for adjusting the play.

Figure 18 shows the adjusting screws with which the access ramp's play can be adjusted. The adjusting screws are provided at both sides of the access ramp.



Figure 18: Adjusting Screws (1)

The inner screws (on the left in Figure 18) are used to adjust the play of the ramp that can be folded down inside the vehicle. The inner screws are only available on an Easy-Flex-Ramp.

If the locking bolt is difficult to close, there is too little play. If the access ramp is loose when in the vertical position or rattles when driving, there is too much play.

To adjust the access ramp's play:

- ⇒ Remove the black protective caps from the adjusting screws.
- ⇒ Loosen the locking nuts (see Figure 19, page 37).



Figure 19: Adjusting Screws (2)

⇒ Turn the screw in to reduce the amount of play.

Or:

⇒ Turn the screw out to increase the amount of play.

⇒ Tighten the locking nut once again.

⇒ Check the access ramp's play once again.

⇒ Repeat the procedure as often as necessary.

⇒ Replace the black protective caps.

8 Inspection Log Book

Before using the access ramp for the first time, it must be inspected by a technical expert.

If the access ramp is used commercially or as a public utility, this inspection must be repeated at intervals of not more than one year

During inspection, faults affecting the safety should be systematically identified and remedial action taken.

The items to be inspected are listed in the inspection list (see Section 8.2, page 40).



A technical expert is someone who, by virtue of their technical training and experience, has sufficient knowledge in the field of access ramps for disabled vehicles and is sufficiently familiar with the relevant national occupational health and safety regulations, accident prevention regulations and recognised rules of sound engineering practice (e.g. BG rules, DIN standards, VDE provisions, technical rules and regulations of other member states of the European Union or Turkey or other signatory states of the Agreement of the European Economic Area) that they are able to assess the safe working condition of such access ramps for disabled vehicles.



8.1 Inspection Log Book Master Data Sheet

Master Data Sheet for Easy-Ramp and Easy-Flex-Ramp

Serial No. /Type*

Owner

.....

.....

Year built:

Commissioned on

*See rating plate

8.2 Inspection List

For a routine, yearly inspection by a technical expert (see information on page 38).

| Mechanical parts | |
|-------------------------|--|
| Fixtures on the vehicle | All securing points in / on the vehicle are in a perfect, tight condition |
| Pivoting elements | All bolts and screwed connections are in a perfect, tight condition |
| General | Functional test and visual inspection of the condition of the entire access ramp |

| Electrical System | |
|-----------------------------|------------------|
| Rocker switch | Functional check |
| Electrically operated belts | Functional check |



8.3 Inspection Results

Results of the inspection prior to commissioning..

To be observed if installed by a third-party company

Installation inspection for

Serial No. /Type*

Installation has been carried out correctly!

Place, date

.....
Installation company / Company stamp

.....
Technical expert / Signature



Inspection result of a yearly inspection

Serial No. /Type*

Sheet No.

The access ramp was subject to a routine inspection on

No / the following defects were found:

.....
.....
.....
.....

There are no issues to preclude continued operation.

Place, date

.....
Company / Company stamp

.....
Signature (technical expert)

All defects have been remedied.

Confirmation by the owner or his representative with date and signature

.....
Place / date

.....
Signature (owner)



Inspection result of a yearly inspection

Serial No. /Type*

Sheet No.

The access ramp was subject to a routine inspection on

No / the following defects were found:

.....
.....
.....
.....

There are no issues to preclude continued operation.

Place, date

.....
Company / Company stamp

.....
Signature (technical expert)

All defects have been remedied.

Confirmation by the owner or his representative with date and signature

.....
Place / date

.....
Signature (owner)



Inspection Log Book

Inspection result of a yearly inspection

Serial No. /Type*

Sheet No.

The access ramp was subject to a routine inspection on

No / the following defects were found:

.....
.....
.....
.....

There are no issues to preclude continued operation.

Place, date

.....
Company / Company stamp

.....
Signature (technical expert)

All defects have been remedied.

Confirmation by the owner or his representative with date and signature

.....
Place / date

.....
Signature (owner)



Inspection result of a yearly inspection

Serial No. /Type*

Sheet No.

The access ramp was subject to a routine inspection on

No / the following defects were found:

.....
.....
.....
.....

There are no issues to preclude continued operation.

Place, date

.....
Company / Company stamp

.....
Signature (technical expert)

All defects have been remedied.

Confirmation by the owner or his representative with date and signature

.....
Place / date

.....
Signature (owner)

9 Electrical Circuit Diagram




Figure 20: Electrical Circuit Diagram

10 Customer Service

The AMF-Brunns customer service department will be more than pleased to assist in ordering spare parts, maintenance and repair work and help with general problems or queries.

The address is:

AMF-Brunns GmbH & Co. KG
Hauptstraße 101
D - 26689 Apen

Tel.: +49 (0) 44 89 / 72 72-30

Fax: +49 (0) 44 89 / 62 45

service.hubmatik@amf-bruns.de

www.amf-bruns.de



NOTE

Guarantee work on the access ramp must only be carried out with the prior agreement of AMF-Brunns GmbH & Co. KG.

The costs of such work will not be accepted by AMF-Brunns without prior agreement.



AMF-Bruns GmbH & Co. KG | Hauptstraße 101 | 26689 Apen
Telephone +49 (0) 4489 / 72 72 30 | Fax +49 (0) 4489 / 62 45
service.hubmatik@amf-bruns.de

www.amf-bruns.de



Reg.-No. Q1 0105027