

Vehicle conversion

PremiumPartnerMAN Truck & Bus

WHEELCHAIR ACCESSIBLE VEHICLE

MAN TGE





www.amf-bruns.de



WHEELCHAIR ACCESSIBLE VEHICLE

MAN TGE



The wheelchair enters the vehicle via lift. AMF-Bruns PROTEKTOR system for wheelchair and occupant restraint.



Smartfloor provides flexibility inside your vehicle.



Optional head and backrest FUTURE**SAFE** for a maximum of safety.

THESE FACTS MAKE ALL THE DIFFERENCE

- » flexible & individual planned Smartfloor Aluminium floor
- » quick & flexible configuration of vehicle's interior with removable seats
- » Smartseats with original upholstery or artificial leather
- » Linearlift AL1 is a fully automatic lift for rear or side door. See-through platform design ensures excellent rear view and meets DIN 75078
- » electric step for easy entry process in the secound and third row
- » Optional head and backrest FutureSafe is reliable, space saving and comfortable. floormounted or wallmounted available
- » interior panelling by AMF-Bruns



Tomorrow's safety today!

FUTURE**SAFE®**

Head- & backrest

FutureSafe is a safe, spacesaving and comfortable patented safety system, which is flexibly customisable and comes with an integrated shoulder belt for the wheelchair occupant.

PROTEKTOR®

Wheelchair & occupant restraint system

The AMF-Bruns Protektor restraint system meets DIN 75078 part II, ISO 10542.

PROSAFE®

"Kraftknotenadapter" Integrated lap belt

With the "Kraftknotenadapter" system the wheelchair is quickly and precisely fitted to the vehicle's floor. It comes with a pelvis belt for the wheelchair occupant.

AMF-BRUNS PROVIDES THE MAXIMUM OF SAFETY

- » The AMF-Bruns 4-point restraint system meets DIN 75078, and ensures that the wheelchairs are safely fixed to the vehicle's floor
- » Linearlift with EMC testing according to EU regulations
- » FutureSafe meets ISO 10542 and DIN 75078-2
- » Vehicle conversion according to DIN 75078-1/2



AMF-Bruns GmbH & Co. KG | Hauptstraße 101 | D-26689 Apen Telefon +49 (0) 44 89 / 72 72 -33 / -50 / -38 | Fax +49 (0) 44 89 / 62 45 hubmatik@amf-bruns.de



