

Safety systems

AMF-BRUNS **PROTEKTOR® 2.0**

Tomorrow's safety today!

Occupant & Wheelchair Restraint System according to ISO 10542





www.amf-bruns.de



AME-BRUNS PROTEKTOR® 2.0

Occupant & Wheelchair Restraint System according to ISO 10542



The lightweight, black ABS cover of the BLACKSERIES retractors is robust and suitable for daily use in your vehicle.

THESE FACTS MAKE ALL THE DIFFERENCE

- ABS cover, durable and lightweight
- · Semi-automatic pull-out
- Rotary knob for additional tensioning without effort
- · Easy, intuitive operation
- · Compact, low design
- Wheelchair attachment with loop, hook or carabiner
- Floor anchorage with 4 stud fitting or mono fitting
- Retractor system includes: 4x retractor belt + 1x pelvis belt (Pro 4) or 2x retractor belt + 2x static belt + 1x pelvis belt (Pro 2)
- Tested safety according to ISO 10542



The noble cover of the SILVERSERIES wheelchair restraints is made of zinc die-casting - a material that makes the retractors particularly robust and which has proven itself in daily use in vehicles.

THESE FACTS MAKE ALL THE DIFFERENCE

- · Zinc die-casting cover, noble and robust
- Fully automatic pull-out
- · Self-tensioning
- · Rotary knob for additional tensioning without effort
- · Easy, intuitive operation
- · Compact, low design
- Wheelchair attachment with loop, hook or carabiner
- Floor anchorage with 4 stud fitting or mono fitting
- Retractor system includes: 4x retractor belt + 1x pelvis belt (Pro 4) or 2x retractor belt + 2x static belt + 1x pelvis belt (Pro 2)
- Tested safety according to ISO 10542

The new generation of the PROTEKTOR® 4-point safety system convinces with its high safety aspect, aesthetic design functionality as well as high-quality materials.

The retractors are anchored to the vehicle floor with a 4 stud fitting or mono fitting and can be easily tethered to the wheelchair with a loop, hook or carabiner attachment.

The wheelchair restraint system is supplemented with a pelvis and shoulder belt to secure the wheelchair occupant.

The BLACKSERIES and SILVERSERIES are compatible with all wheelchairs and vehicles and is ideal for use of wheelchairs with an integrated lap belt PROSAFE.



Tomorrow's safety today!

FUTURE**SAFE®**

Head & Backrest

FUTURE**SAFE** is a safe, space-saving and com-fortable patented safety system, which is flexibly customisable and comes with an integrated shoulder belt for the wheelchair occupant.

PROTEKTOR®

Occupant & Wheelchair **Restraint System**

The PROTEKTOR® restraint system from AMF-Bruns meets DIN 75078 part II, ISO 10542.

PROSAFE®

"Kraftknotenadapter" Integrated Lap Belt

With the "Kraftknotenadapter" system the wheelchair is quickly and precisely fitted to the vehicle's floor. It comes with a pelvic belt for the wheelchair occupant.



AMF-Bruns GmbH & Co. KG | Hauptstraße 101 | D-26689 Apen Phone +49 (0) 44 89 / 72 72 50 | Fax +49 (0) 44 89 / 62 45 hubmatik@amf-bruns.de

www.amf-bruns.de

