



Lift systemes

# AMF-BRUNS CASSETTE LIFT K70

Premium**Partner**  
of Volkswagen Commercial Vehicles

Volkswagen Multivan T7



Tomorrow's  
**safety** today!



[www.amf-bruns.de](http://www.amf-bruns.de)



# AMF-BRUNS CASSETTE LIFT K70

The ultra-slim AMF-Bruns cassette lift for the side sliding door



The lateral access with the flat side edge of the lifting platform saves space when driving up with the wheelchair.



The ultra-slim K70 cassette lift combined with the vehicle heightening offers the best possible ground clearance.

## THESE FACTS MAKE ALL THE DIFFERENCE

- » Load capacity 300 kg, complies with EN 1756, II
- » Ultra-slim design (70 mm) guarantees ultimate ground clearance
- » Illuminated controls, waterproof and shockproof
- » Automatic roll-off barrier (also operates on uneven roads and snow)
- » Powerful, silent running function ensures smooth operation of the lifting platform
- » Large lifting platform with side access; also suitable for large electric wheelchairs
- » Fully closed stainless steel cassette offers safe protection of the lift system, even in extreme weather conditions
- » No changes to the car body necessary
- » Vehicle interior remains unchanged
- » Vehicle heightened by approx. 55 mm without impairing driving comfort
- » Manual emergency operation with hand pump
- » EMC tested according to EU regulations

## TECHNICAL DATA

» Height	70 mm	» Deploying/ stowing	approx. 12 sec.
» Load capacity	300 kg	» Control	By push buttons
» Platform size	approx. 760 x 1,200 mm	» Power	Lifting and lowering electro-hydraulic 12 V; deploying and stowing by rack and pinion drive 12 V
» Unladen weight	approx. 170 kg	» Optional	Bluetooth remote control including app operation (iOS and Android)
» Lifting time	approx. 7 sec.		
» Lowering time	approx. 7 sec.		

## MAXIMUM SAFETY

- » Automatic roll-over flap prevents feet from being trapped when the platform is lifted
- » Fully automatic roll-over protection barrier
- » Two LED flashing lights that switch on automatically to warn in poor visibility conditions
- » Safety switches prevent incorrect operation while driving



AMF-Bruns GmbH & Co. KG | Hauptstraße 101 | D-26689 Apen  
 Telefon +49 (0) 44 89 / 72 72 50 | Fax +49 (0) 44 89 / 62 45  
 hubmatik@amf-bruns.de

[www.amf-bruns.de](http://www.amf-bruns.de)



DIN EN ISO 9001  
 REG.-NR. 01 0105027